

# Kartell

Labware division

## dispolab®

### Micro-pipette Tip Compatibility Charts



#### **R & L Slaughter Ltd.**

Unit 14, Carnival Park, Carnival Close, Basildon, Essex. SS14 3WN  
Tel: 01708 227140 Fax: 01708 228728 email: [sales@slaughter.co.uk](mailto:sales@slaughter.co.uk)  
[www.slaughter.co.uk](http://www.slaughter.co.uk)

# Kartell - dispolab®

## Micro-Pipette Tip Compatibility Chart - Brand Order

Pipettor Brand	Capacity	Colour	Loose [Bulk] (Cat. No.)	Boxed Racks (Cat. No.)	Refill Racks (Cat. No.)	Boxed Racks Economic (Cat. No.)	Stacked Racks (Cat. No.)	Sterile Filter Tip (Cat. No.)
Beckman®	- 1000 µl	Neutral	965	--	--	--	--	--
Biohit®	0.5 - 10 µl	Crystal	956	9646	9747	--	--	9638
Biohit®	2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632 100 µl neutral
Biohit®	2 - 300 µl	Neutral	971	9644	9752	--	--	--
Biohit®	10 - 300 µl	Neutral	935 *	9647 *	9754 *	--	--	9637
Biohit®	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Biohit®	200 - 1000 µl	Blue	957	--	--	9615	9628	--
Biohit®	1000 - 5000 µl	Neutral	964	9649	--	--	--	--
Brand®	2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632 100 µl neutral
Brand®	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Brand®	1000 - 5500 µl	Neutral	942	9847	--	--	--	--
Eppendorf®	0.5 - 10 µl	Crystal	959	9650	9746	9624	--	9630
Eppendorf®	2 - 20 µl	Neutral	979	9657	9745	--	--	9639
Eppendorf®	2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632 100 µl neutral
Eppendorf®	10 - 300 µl	Neutral	935 *	9647 *	9754 *	--	--	9637
Eppendorf®	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Eppendorf®	200 - 1000 µl	Blue	957	--	--	9615	9628	--
Eppendorf®	1000 - 5000 µl	Neutral	964	9649	--	--	--	--
Finnpipette®	0.5 - 10 µl	Crystal	956	9646	9747	--	--	9638
Finnpipette®	2 - 300 µl	Neutral	971	9644	9752	--	--	--
Finnpipette®	1000 - 5500 µl	Neutral	942	9847	--	--	--	--
Finnpipette®	1000 - 10000 µl	Neutral	943	9648	--	--	--	--
Gilson® (P-2 + P-10)	0.5 - 10 µl	Crystal	956	9646	9747	--	--	9638
Gilson® (P-20)	0.5 - 10 µl	Crystal	959	9650	9746	9624	--	9630
Gilson® (P-20 to 200)	2 - 200 µl	Yellow	958	9642	9750	9585	9627	--
Gilson® (P-20 to 200)	5 - 200 µl	Neutral	940 *	9643 *	9751 *	--	--	9635 100 µl 9636 160 µl
Gilson® (P-500 + 1000)	100 - 1000 µl	Blue	941	9645	9753	--	--	9633 neutral
Gilson® (Universal)	2 - 20 µl	Neutral	979	9657	9745	--	--	--
Gilson® (Universal)	2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632 100 µl neutral
Gilson® (Universal)	2 - 300 µl	Neutral	971	9644	9752	--	--	--
Gilson® (Universal)	10 - 300 µl	Neutral	935 *	9647 *	9754 *	--	--	9637
Gilson® (Universal)	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Gilson® (Universal)	200 - 1000 µl	Blue	957	--	--	9615	9628	--
Gilson® (Universal)	1000 - 5000 µl	Blue Tinge	984	9648	--	--	--	--

Continues on top of the next page



# Kartell - dispolab®

## Micro-Pipette Tip Compatibility Chart - Brand Order Cont.

Pipettor Brand	Capacity	Colour	Loose [Bulk]	Boxed Racks	Refill Racks	Boxed Racks Economic	Stacked Racks	Sterile Filter Tip
			(Cat. No.)	(Cat. No.)	(Cat. No.)	(Cat. No.)	(Cat. No.)	(Cat. No.)
Gilson® (Universal)	1000 - 10000 µl	Neutral	943	9848	--	--	--	--
Kartell®	0.5 - 10 µl	Crystal	959	9650	9746	9624	--	9630
Kartell®	2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632 100 µl neutral
Kartell®	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Kartell®	200 - 1000 µl	Blue	957	--	--	9615	9628	--
Kartell®	1000 - 5000 µl	Neutral	964	9649	--	--	--	--
Kartell®	1000 - 10000 µl	Neutral	943	9848	--	--	--	--
MLA®	5 - 200 µl	Neutral	968	--	--	--	--	--
MLA®	200 - 1000 µl	Neutral	969	--	--	--	--	--
Nichiryo®	0.5 - 10 µl	Crystal	959	9650	9746	9624	--	9630
Nichiryo®	2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632 100 µl neutral
Nichiryo®	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Nichiryo®	1000 - 5000 µl	Blue Tinge	984	9648	--	--	--	--
Oxford Slimline®	5 - 200 µl	Neutral	966	--	--	--	--	--
Oxford Slimline®	250 - 1000 µl	Blue	967	--	--	--	--	--
Oxford® (orig.)	5 - 200 µl	Neutral	962	--	--	--	--	--
Oxford® (orig.)	250 - 1000 µl	Green	963	--	--	--	--	--
Oxford® (orig.)	1000 - 5000 µl	Neutral	964	9649	--	--	--	--
Sherwood® / Lancer®	- 1000 µl	Neutral	965	--	--	--	--	--
Socorex®	0.5 - 10 µl	Crystal	959	9650	9746	9624	--	9630
Socorex®	2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632 100 µl neutral
Socorex®	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Socorex®	1000 - 5000 µl	Blue Tinge	984	9648	--	--	--	--
Socorex®	1000 - 10000 µl	Neutral	943	9848	--	--	--	--
Treff®	100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631 neutral
Treff®	1000 - 5500 µl	Neutral	942	9847	--	--	--	--

### Legend:

\* suffix = Tip is graduated

Filter Tips are 'Rack Packed' AND Gamma Ray Sterilised

for Labsystem® refer Finnpiptette®

for Pipettman® refer Gilson®

## Micro-Pipette Tip Compatibility Chart - Capacity Order

Capacity	Colour	Loose [Bulk] (Cat No.)	Boxed Racks (Cat. No.)	Refill Racks (Cat. No.)	Boxed Racks Economic (Cat. No.)	Stacked Racks (Cat. No.)	Sterile Filter Tip (Cat. No.)	Compatibility
0.5 - 10 µl	Crystal	956	9646	9747	--	--	9638	Biohit <sup>®</sup> , Finnpiquette <sup>®</sup> , Gilson <sup>®</sup> (P-2 + P-10)
0.5 - 10 µl	Crystal	959	9650	9746	9624	--	9630	Eppendorf <sup>®</sup> , Gilson <sup>®</sup> (P-20), Kartell <sup>®</sup> Pluripet PL 10, Nichiryo <sup>®</sup> , Socorex <sup>®</sup>
2 - 20 µl	Neutral	979	9657	9745	--	--	9639 (Eppendorf <sup>®</sup> )	Eppendorf <sup>®</sup> , Gilson <sup>®</sup>
2 - 200 µl	Yellow	960 980 *	9641	9749	9605	9626	9632	Biohit <sup>®</sup> , Brand <sup>®</sup> , Eppendorf <sup>®</sup> , Gilson <sup>®</sup> , Kartell <sup>®</sup> , Nichiryo <sup>®</sup> , Socorex <sup>®</sup> , Kartell Multichannel
2 - 200 µl	Yellow	958	9642	9750	9585	9627	--	Gilson <sup>®</sup> (P-20 to P-200), Kartell Multichannel
5 - 200 µl	Neutral	940 *	9643 *	9751 *	--	--	9635 9636	Gilson <sup>®</sup> (P-20 to P-200), Kartell Multichannel
5 - 200 µl	Neutral	968	--	--	--	--	--	MLA <sup>®</sup>
5 - 200 µl	Neutral	966	--	--	--	--	--	Oxford <sup>®</sup> Slimline
2 - 300 µl	Neutral	971	9644	9752	--	--	--	Biohit <sup>®</sup> , Gilson <sup>®</sup> , Finnpiquette <sup>®</sup>
10 - 300 µl	Neutral	935 *	9647 *	9754 *	--	--	9637	Eppendorf <sup>®</sup> , Gilson <sup>®</sup> , Biohit <sup>®</sup>
100 - 1000 µl	Blue	961 981 *	9640	9748	--	--	9631	Biohit <sup>®</sup> , Brand <sup>®</sup> , Gilson <sup>®</sup> , Kartell <sup>®</sup> , Nichiryo <sup>®</sup> , Socorex <sup>®</sup>
100 - 1000 µl	Blue	941	9645	9753	--	--	9633	Gilson <sup>®</sup> (P-500 + P-1000)
200 - 1000 µl	Blue	957	--	--	9615	9628	--	Biohit <sup>®</sup> , Eppendorf <sup>®</sup> , Gilson <sup>®</sup> , Kartell <sup>®</sup>
200 - 1000 µl	Neutral	969	--	--	--	--	--	MLA <sup>®</sup>
250 - 1000 µl	Green	963	--	--	--	--	--	Oxford <sup>®</sup> (orig.)
250 - 1000 µl	Blue	967	--	--	--	--	--	Oxford <sup>®</sup> Slimline
- 1000 µl	Neutral	965	--	--	--	--	--	Beckman <sup>®</sup> , Sherwood <sup>®</sup> / Lancer <sup>®</sup>
1000 - 5000 µl	Blue Tinge	984	9648	--	--	--	--	Gilson <sup>®</sup> , Nichiryo <sup>®</sup> , Socorex <sup>®</sup> , Kartell <sup>®</sup> Pluripet std & Pluripet II
1000 - 5000 µl	Neutral	964	9649	--	--	--	--	Kartell <sup>®</sup> , Eppendorf <sup>®</sup> , Biohit <sup>®</sup> , Oxford <sup>®</sup> , Kartell <sup>®</sup> Pluripet std
1000 - 5500 µl	Neutral	942	9847	--	--	--	--	Finnpiquette <sup>®</sup> , Treff <sup>®</sup> , Brand <sup>®</sup>
1000 - 10000 µl	Neutral	943	9848	--	--	--	--	Gilson <sup>®</sup> , Kartell <sup>®</sup> , Socorex <sup>®</sup> , Finnpiquette <sup>®</sup> , Gilson <sup>®</sup>

### Legend:

\* suffix = Tip is graduated  
 Filter Tips are 'Rack Packed' AND Gamma Ray Sterilised  
 for Labsystem<sup>®</sup> refer Finnpiquette<sup>®</sup>  
 for Pipettman<sup>®</sup> refer Gilson<sup>®</sup>