

Sharps containers, Code YELLOW, SHARPAK™



All SHARPAK™ Code YELLOW sharps containers feature yellow coloured closures and are recommended for the segregation, collection and disposal of partially discharged sharps waste which are NOT contaminated with cytotoxic or cytostatic medicines. Code YELLOW SHARPAK™ sharps containers MUST only be disposed of by incineration. Tough and durable. Cost effective.

All SHARPAK™ Code YELLOW sharps containers are approved to BS7320:1990 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified.

Please note: For disposal of sharps contaminated with cytotoxic or cytostatic medicines please refer to sharps containers, code PURPLE, SHARPAK™.

- Yellow colour coded closures to meet the recommendations on waste segregation as detailed in the Department of Health's technical memorandum HTM-07-01
- Secure closure with locking facility
- Distinctive colouring and clear labelling
- Sequential numbering for traceability
- Supplied assembled

Capacity (l)	WxDxH (mm)	Type	Pk	Cat. No.
3,5	177x206x195	Sharpak™ 35 PLUS, yellow lid	1	129-1060
6	220x234x229	Sharpak™ 60 PLUS, yellow lid	1	129-1058
12	258x268x325	Sharpak™ 120 PLUS, yellow lid	1	129-1052
22	315x327x383	Sharpak™ 220 PLUS, yellow lid	1	129-1055

Sharps containers, Code ORANGE, SHARPAK™



All SHARPAK™ Code ORANGE sharps containers feature orange coloured closures and are recommended for the segregation, collection and disposal of used and fully discharged sharps waste which is NOT contaminated with cytotoxic or cytostatic medicines or sharps NOT contaminated with medicinal products. Code ORANGE SHARPAK™ sharps containers can be disposed of by incineration or other approved method. Tough and durable. Cost effective.

All SHARPAK™ Code ORANGE sharps containers are approved to BS7320:1990 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified.

Please note: For disposal of sharps contaminated with cytotoxic or cytostatic medicines please refer to SHARPAK™, Code PURPLE containers.

- Orange colour coded closures to meet the recommendations on waste segregation as detailed in the Department of Health's technical memorandum HTM-07-01
- Secure closure with locking facility
- Distinctive colouring and clear labelling
- Sequential numbering for traceability
- Supplied assembled

Capacity (l)	WxDxH (mm)	Type	Pk	Cat. No.
3,5	177x206x195	Sharpak 35 PLUS, orange lid	1	129-1059
6	220x234x229	Sharpak 60 PLUS, orange lid	1	129-1057
12	258x268x325	Sharpak 120 PLUS, orange lid	1	129-1051
22	315x327x383	Sharpak 220 PLUS, orange lid	1	129-1054

ALL YOU NEED FOR CHEMICAL ANALYSIS

A great companion for chemists in any industry

- Sampling and sample preparation
- Spectroscopy, chromatography and physical measurement
- Waste and hazard management

Request your copy from your local VWR sales office or vwr.com