



Magnetic stirring bar retrievers



PTFE covered, strong Alnico V magnetic core. For the removal of stirring magnets from vessels of all kinds.

- Very good chemical resistance
- Temperature resistant, suitable for use at high and low temperatures (–200 to +280 °C)

Length (mm)	Ø (mm)	Pk	Cat. No.
150	10	1	442-0523
250	10	1	442-0524
350	10	1	442-0525
450	10	1	442-0526



Magnetic stirring bar retrievers



PTFE covered, rare earth magnetic core. For the removal of stirring magnets from vessels of all kinds.

- Very good chemical resistance
- Temperature resistant, suitable for use at high and low temperatures (–200 to +280 °C)

Length (mm)	Ø (mm)	Pk	Cat. No.
150	10	1	442-4549
250	10	1	442-4550
350	10	1	442-4551
450	10	1	442-4552

ALL YOU NEED FOR BIOBANKING

- Sample storage
- Sample handling
- Safety

Request your copy from your local VWR sales office or vwr.com