

Pyridine GPR RECTAPUR®

Assay.....	Min. 99 %
Evaporation residue	Max. 20 ppm
Water	Max. 0.5 %

Cat. No.	Pk	Pack type
8226.2500	2,5 l	Glass bottle

Pyridine TECHNICAL

Assay.....	Min. 99 %
------------	-----------

Cat. No.	Pk	Pack type
27194.295	1 l	Glass bottle

Pyrogallic acid

See Pyrogallol..... p.400

Pyrogallol

1,2,3-Trihydroxybenzene , Pyrogallic acid

Warning

H341 H302+H312+H332 H412
P201 P281 P273 P302+P352 P304+P340 P309+P311

CAS 87-66-1

Index 604-009-00-6

EINECS: 201-762-9

UN: 2811

ADR 6.1,III

C₆H₅O₃

M.W. 126.11 g/mol

Density: 1.453 g/cm³ (20 °C)

Boiling Pt: 309 °C (1013 hPa)

Melting Pt: 130 to 132 °C

Storage Temperature: Ambient temperature



Pyrogallol AnalR NORMAPUR® analytical reagent

Assay.....	Min. 99.5 %	Acidity.....	Max. 0.005 meq/g
Melting point.....	131 to 134 °C	Heavy metals (as Pb)	Max. 5 ppm
Ignition residue (SO ₄).....	Max. 50 ppm	Cl (Chloride)	Max. 10 ppm
SO ₄ (Sulphate).....	Max. 50 ppm	Fe (Iron).....	Max. 10 ppm

Cat. No.	Pk	Pack type
27206.231	250 g	Glass bottle for solids

Pyrogallol GPR RECTAPUR®

IR Spectrum.....	Passes test
Melting point.....	131 to 133 °C
Ignition residue (SO ₄).....	Max. 0.1 %
Water	Max. 1.0 %

Cat. No.	Pk	Pack type
20647.237	250 g	Glass bottle for solids

Pyronine G

3,6-Bis(dimethylamino)xanthylium chloride , Pyronine Y

CAS 92-32-0

EINECS: 202-147-8

C₁₇H₁₉ClN₂O

M.W. 302.8 g/mol

Melting Pt: 250 to 260 °C

Storage Temperature: Ambient temperature

VWR CHEMICALS BDH PYRONINE G, ultrapure

Used in Unna-Pappenheim stain to demonstrate RNA (red) in contrast to DNA (green). Also used for bacteria staining.

Lambda max (0.1 g/500 ml then 1:100, Water)..... 543 - 553 nm

Cat. No.	Pk	Pack type
0207-10G	10 g	Glass bottle

Pyruvic acid sodium salt

See Sodium pyruvate..... p.455



VWR BDH PROLABO® CHEMICALS

HIPERSOLV® CHROMANORM®

- High purity solvents for HPLC applications
- Designed to meet your requirements in analysis and quality control



VWR BDH PROLABO® CHEMICALS

BIOCHEMICALS, MORE THAN 100 YEARS EXPERIENCE

Please contact your local VWR sales office for more information.



VWR BDH PROLABO® CHEMICALS

REAGENTS FOR HISTOPATHOLOGY

Complete range of reagents for cell diagnostics including the new Q Path® range