




EDTA TUBES

PROVETTE CON EDTA






K2 EDTA TUBES PROVETTE CON K2 EDTA

STERILE R

CODE codice	SIZE dimensioni	DESCRIPTION descrizione	DRAWING aspirazione	COLOUR colore		SHELF-LIFE scadenza	PACKAGING confezion.
13501	13 x 75 mm	K2 EDTA K2 EDTA	1 ml	Lavender Viola		18 months mesi	100 / 1000
13505	13 x 75 mm	K2 EDTA K2 EDTA	2 ml	Lavender Viola		18 months mesi	100 / 1000
13520	13 x 75 mm	K2 EDTA K2 EDTA	2 ml	Translucent Lavender Viola trasparente		18 months mesi	100 / 1000
135209	13 x 75 mm	K2 EDTA K2 EDTA	2 ml	Black Nero		18 months mesi	100 / 1000
13510	13 x 75 mm	K2 EDTA K2 EDTA	3 ml	Lavender Viola		18 months mesi	100 / 1000
135100	13 x 75 mm	K2 EDTA K2 EDTA	3 ml	Translucent Lavender Viola trasparente		18 months mesi	100 / 1000
135201	13 x 75 mm	K2 EDTA K2 EDTA	3 ml	Black Nero		18 months mesi	100 / 1000
135300	13 x 75 mm	K2 EDTA K2 EDTA	4 ml	Lavender Viola		18 months mesi	100 / 1000
135306	13 x 75 mm	K2 EDTA K2 EDTA	4 ml	Translucent Lavender Viola trasparente		18 months mesi	100 / 1000
135171	13 x 75 mm	K2 EDTA K2 EDTA	4 ml	Pink Rosa		18 months mesi	100 / 1000
135400	13 x 100 mm	K2 EDTA K2 EDTA	6 ml	Lavender Viola		18 months mesi	100 / 1000
135420	13 x 100 mm crossmatch	K2 EDTA K2 EDTA	6 ml	Fuchsia Fucsia		18 months mesi	100 / 1000
135560	13 x 100 mm	K2 EDTA K2 EDTA	6 ml	Light green Verde acqua		18 months mesi	100 / 1000
13560	16 x 100 mm	K2 EDTA K2 EDTA	9 ml	Lavender Viola		18 months mesi	100 / 1000


K3 EDTA TUBES
PROVETTE CON K3 EDTA

CODE <i>codice</i>	SIZE <i>dimensioni</i>	DESCRIPTION <i>descrizione</i>	DRAWING <i>aspirazione</i>	COLOUR <i>colore</i>	SHELF-LIFE <i>scadenza</i>	PACKAGING <i>confezion.</i>	
13001	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	1 ml	Lavender <i>Viola</i>		18 months <i>mesi</i>	100 / 1000
13005	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	2 ml	Lavender <i>Viola</i>		18 months <i>mesi</i>	100 / 1000
13840	13 x 75 mm IN	K3 EDTA <i>K3 EDTA</i>	2 ml	Translucent Lavender <i>Viola trasparente</i>		18 months <i>mesi</i>	100 / 1000
130200	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	2 ml	Translucent Lavender <i>Viola trasparente</i>		18 months <i>mesi</i>	100 / 1000
13148	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	2 ml	Black <i>Nero</i>		18 months <i>mesi</i>	100 / 1000
13010	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Lavender <i>Viola</i>		18 months <i>mesi</i>	100 / 1000
13125	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Translucent Lavender <i>Viola trasparente</i>		18 months <i>mesi</i>	100 / 1000
13150	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Black <i>Nero</i>		18 months <i>mesi</i>	100 / 1000
13161	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	White <i>Bianco</i>		18 months <i>mesi</i>	100 / 1000
131100	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Yellow <i>Giallo</i>		18 months <i>mesi</i>	100 / 1000
13025	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Dark green <i>Verde scuro</i>		18 months <i>mesi</i>	100 / 1000
13136	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Orange <i>Arancio</i>		18 months <i>mesi</i>	100 / 1000
13131	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Fuchsia <i>Fucsia</i>		18 months <i>mesi</i>	100 / 1000
13182	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	3 ml	Rusty <i>Terra di Siena</i>		18 months <i>mesi</i>	100 / 1000
13030	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	4 ml	Lavender <i>Viola</i>		18 months <i>mesi</i>	100 / 1000
13220	13 x 75 mm	K3 EDTA <i>K3 EDTA</i>	4 ml	Dark blue <i>Blu</i>		18 months <i>mesi</i>	100 / 1000
13040	13 x 100 mm	K3 EDTA <i>K3 EDTA</i>	6 ml	Lavender <i>Viola</i>		18 months <i>mesi</i>	100 / 1000
13060	16 x 100 mm	K3 EDTA <i>K3 EDTA</i>	9 ml	Lavender <i>Viola</i>		18 months <i>mesi</i>	100 / 1000

Recommended use**Mixing indications:** immediately after blood collection, gently invert the sample 6-8 times.**Centrifugation:** do not centrifuge.**Sample preservation:**

Storage temperatures	Up to 24 °C	at 2-4 °C
Maximum preservation time	6 hours	≤ 24 hours

Before processing the haematological test mix carefully**To obtain plasma:** centrifuge at 1300 g for 10 minutes at 20-25°C.**Sample preservation:** for plasma separated from blood cells see following table:

Storage temperatures	Up to 24 °C	at 2-4 °C	≤ -20 °C
Maximum preservation time	≤ 8 hours	> 8 hours ≤ 48 hours	> 48 hours

Raccomandazioni d'uso**Indicazioni per miscelare:** subito dopo il prelievo agitare 6-8 volte il campione per inversione lenta.**Centrifugazione:** non centrifugare.**Conservazione del campione:**

Temperature di conservazione	Fino a 24°C	a 2-4 °C
Tempo massimo di conservazione	6 ore	≤ 24 ore

Prima di effettuare test ematologici miscelare accuratamente**Per ottenere plasma:** centrifugare a 1300 g per 10 minuti a 20-25°C.**Conservazione del campione:** seguire le indicazioni della tabella sottostante per il plasma separato dalla parte corpuscolata:

Temperature di conservazione	Fino a 24°C	a 2-4 °C	≤ -20 °C
Tempo massimo di conservazione	≤ 8 ore	> 8 ore ≤ 48 ore	> 48 ore