

## Cell culture tubes, Nunclon™∆ Thermo Scientific



### PS with PE caps, sterile

Round bottom tubes with screw or push-on cap for standard cell culture, or flat-sided tubes for the culture of adherent cells.

- Flat side allows easy microscopy a coverslip can be used in the tube
- Medium retained in the flat-sided tube in the horizontal position
- Excellent optical quality
- Certified surface treatment for optimal cell attachment and growth

Capacity (ml)	Description	Pk	Cat. No.
7	Round tubes with push-on cap	600	734-2036
7	Round tubes with screw cap	600	734-2037
3	Flat-sided tubes with screw cap	450	734-2068

## Sample tubes, disposable, Nunc™ Thermo Scientific



#### PS

PS tubes are suitable for a broad range of applications.

- Graduation marked
- PS tubes are transparent
- Stoppers sold separately

Length (mm)	Ø ext. (mm)	Description	Material	Sterile	Pk	Cat. No.
70	11	Round bottom	PS	-	3.600	734-0493
110	17	Conical bottom	PS	-	1.800	525-0060
					Pk	Cat. No.
Stoppers, PE, for tubes 11 mm external diameter						391-8056
Stoppers, PE, for tubes 12 mm external diameter						732-2680
Stoppers, PE, for tubes 15 mm external diameter						
	70 110 bes 11 mm external diame ses 12 mm external diame	70 11 110 17  bes 11 mm external diameter ses 12 mm external diameter	70 11 Round bottom 110 17 Conical bottom  bes 11 mm external diameter es 12 mm external diameter	70 11 Round bottom PS 110 17 Conical bottom PS  bes 11 mm external diameter  ses 12 mm external diameter	70 11 Round bottom PS - 110 17 Conical bottom PS -  wes 11 mm external diameter  wes 12 mm external diameter	70         11         Round bottom         PS         -         3.600           110         17         Conical bottom         PS         -         1.800           Pk           sees 11 mm external diameter         3.600           sees 12 mm external diameter         3.000



# **ALL YOU NEED FOR CELL CULTURE**

Serum, Media, Reagents, Transfection, RNAi, Electroporation...

Request your copy from your local VWR sales office or vwr.com

