

Continued from previous page

Description	Colour	Well volume (µl)	No. of wells	Sterile	Pk	Cat. No.
<b>Standard Stripwell microplates</b>						
Standard microplate, medium binding*	Clear	360	96	-	100	734-1515
Standard microplate, high binding*	Clear	360	96	-	100	734-1514
Standard microplate, high binding**	Clear	360	96	-	800	732-5504
Standard microplate, TC-treated, sterile***	Clear	360	96	+	50	732-5569
Standard microplate, medium binding	White	360	96	-	100	732-5545
Standard microplate, high binding	Black	360	96	-	100	732-5548
Standard microplate, medium binding	Black	360	96	-	100	732-5546
Standard microplate, high binding	White	360	96	-	100	732-5547

\*Product has a certified surface chemistry  
\*\*Individual 1x8 strips without frame, bulk packed  
\*\*\*Microplates individually packaged with lid

Nunc-Immuno™ stick  
Thermo Scientific



PP tube and stick, PS paddle

Developed for ELISA and RIA tests in the field.

- Gives reliable, semi-quantitative results with built-in negative control
- Hydrophilic MaxiSorp™ paddle surface
- PP tube for minimal adsorption of assay components

Paddle surface area: 5,2 cm²

External LxW (including lid): 49x12 mm

Description	Pk	Cat. No.
Nunc-Immuno™ stick with MaxiSorp™ paddle surface	1.800	391-8283

Tubes, Nunc-Immuno™, Nunc™  
Thermo Scientific



PS with a choice of MaxiSorp™ (hydrophilic) or PolySorp™ (hydrophobic) surfaces

Nunc-Immuno™ tubes are used for solid phase immunoassay methods such as IRMA, ELISA and ILMA methods.

- Nunc StarTubes™ with an increased internal surface, give higher sensitivity and faster assay times
- Very good optical quality and MaxiSorp™ versions have certified binding homogeneity
- Compatible with standard equipment

Capacity (ml)	Description	Ø ext.xH (mm)	Pk	Cat. No.
4,0	Tubes MaxiSorp™	11x70	1.800	735-0019
5,0	Tubes MaxiSorp™	12x75	3.000	735-0018
5,0	Nunc StarTubes™, MaxiSorp™	12x75	3.000	391-8281
5,0	Nunc StarTubes™, PolySorp™	12x75	3.000	735-0072



VWR AND AMRESKO DELIVER  
A NEW LIFE SCIENCE BRAND

Same Pure Product | New Pure Brand

AMRESKO products may not be available in every country,  
please contact your local VWR sales office