


Separation and Concentration

Centrifuges and Sedimentation

Centrifuge Tubes

Conical-Bottom Centrifuge Bottle Adapters

Nalgene



PC, autoclavable

- Nalgene® conical-bottom centrifuge bottles (175 ml) must be used with this adapter
- Fits most commonly used low- and high-speed rotors
- Not interchangeable with other adapters



Ø (mm)	for (ml)	Pk	Cat. No.
61	175	4	525-2400

Round-Bottom Centrifuge Bottle Adapters

Nalgene



PC, autoclavable

- For use in standard centrifuge rotors

Total height of Nalgene® centrifuge bottle with closure and adapter: 137.5 mm



Ø (mm)	for (ml)	Pk	Cat. No.
61.7	250	4	525-2401

Centrifuge Bottle Adapters

Nalgene



LDPE

- Must be used for 250 ml centrifuge bottles with Sorvall GSA rotors
- Facilitates the removal of the bottle from the rotor after centrifugation

Ø (mm)	for (ml)	Pk	Cat. No.
61.5	250	6	525-2300

Centrifuge Tubes Closure

Nalgene



PP, autoclavable

- Designed for both Nalgene® round bottom and conical centrifuge tubes
- A tab facilitates the handling of the closure
- May not provide liquid-tight seal during high-speed centrifugation

Ordering Information: 525-2937 is for 50 ml tubes with rim and round bottom, 525-2938 is for 50 ml round bottom, rimless tubes.



Ø (mm)	for (ml)	Pk	Cat. No.
16	12	20	525-2932
17	15	20	525-2933
18	16	20	525-2934
25	38	20	525-2936
29	50	20	525-2937
30	50	20	525-2938
32	100	20	525-2939