

Special Applications

Healthcare

Haematology

Haematocrit Sealing Compound

Brand



Non-drying vinyl putty on a plastic tray. For rapid sealing of several hundred micro haematocrit capillaries. Up to 24 capillaries can be placed vertically in numbered positions along the sides of the tray.



- Non-drying putty

Description	Pk	Cat. No.
Sealing compound	10	521-9105

Reader, Micro-Haematocrit



For determination of packed cell volume, using 75 mm capillary tubes with total fluid column heights between 40 and 70 mm.

Description	Pk	Cat. No.
Reader	1	720-0100

Haemocytometer, Improved Neubauer

Haemocytometer, double cell, but with Improved Neubauer ruling, to BS 748.

Description	Pk	Cat. No.
Haemocytometer	1	720-0104

Cristaseal



For sealing capillary tubes. The dry end of the filled tube is pressed into the Cristaseal with a twisting motion producing a firm seal with flat boundaries. The tray may be used as a carrier for tubes prior to centrifuging.

- Quicker than flame sealing which can haemolyse the blood

Description	Pk	Cat. No.
Cristaseal	10	720-0101