

Scandium standard solution, 1,000 mg/l Sc in dil. nitric acid

Warning
H319 H315
P280 P302+P352 P305+P351+P338

CAS 7440-20-2
EINECS: 231-129-2
UN: 3264
ADR 8,III

Sc
M.W. 44.96 g/mol
Storage Temperature: Ambient temperature



Scandium standard solution, 1,000 mg/l Sc in dil. nitric acid (from Sc) ARISTAR® standard for ICP-MS

Sc in 2 % HNO₃

Traceable to SRM from NIST, tested in ISO Guide 34 /ISO17025 accredited laboratory.

Cat. No.	Pk	Pack type
456872V	100 ml	Plastic bottle

Supplied with certificate of analysis.

Scandium standard solution, 1,000 mg/l Sc in dil. nitric acid (from Sc₂O₃) ARISTAR® standard for ICP

Sc₂O₃ in HNO₃

Traceable to SRM from NIST, tested in ISO Guide 34 /ISO17025 accredited laboratory.

Cat. No.	Pk	Pack type
455962S	100 ml	Plastic bottle
455964U	500 ml	Plastic bottle

Supplied with certificate of analysis.

NEW Scandium standard solution, 1,000 mg/l Sc in dil. nitric acid AVS TITRINORM® standard for AAS

Cat. No.	Pk	Pack type
86708.180	100 ml	Plastic bottle
86708.260	500 ml	Plastic bottle

Supplied with certificate of analysis.

Schiff's reagent (Feulgen stain)

Density: 1.1 g/cm³ (20 °C)


Schiff's reagent (Feulgen stain) for analysis of aldehydes

Remark: In the particular case of formaldehyde, use an ethyl alcohol solution so that the final solution contains about 50% of ethyl alcohol.

Density (20/4)..... 1.005 to 1.015

Cat. No.	Pk	Pack type
30969.261	500 ml	Glass bottle

Schiff's reagent (Feulgen stain) for staining acc. to Feulgen



IVD

Identity (colour reaction)..... Passes test
Suitability for microscopy..... Passes test
Suitability for electrophoresis..... Passes test

Cat. No.	Pk	Pack type
351204L	500 ml	Glass bottle

Technical data sheet and instructions available on vwr.com

Schiff's staining solution

Storage Temperature: Ambient temperature

NEW Schiff's staining solution Q PATH® for microscopy



Ready to use Schiff's solution for histological staining in easy to use and safe PET bottle.

IVD

Cat. No.	Pk	Pack type
10047025.	6 x 450 ml	Plastic bottle

IVD registered. Instructions for use on vwr.com - just search for the product.

Shorr staining solution

Danger
H225
P210 P243 P280
UN: 1170
ADR 3,II

Storage Temperature: Ambient temperature

