



## Cell counters, Bürker Marienfeld



### Glass

These cell counters for investigations of cell suspensions are equipped with two sets of grids. They are engraved into two ground and polished surfaces on the central ridge (=chamber bottom). This central ridge is located between two elevated ridges that are also ground and polished.

- *In vitro* diagnostic (IVD) applications

Conform to DIN 12847 and the German Calibration Ordinance

**Delivery information:** Supplied in a transparent plastic box.

Description	Pk	Cat. No.
Cell counting chamber, bright lined, Bürker	1	631-0921
Cell counting chamber, bright lined, Bürker-Türk	1	630-1544
Cell counting chamber, bright lined, with clamps, Bürker	1	630-1542
Cell counting chamber, bright lined, with clamps, Bürker-Türk	1	630-1545
Cell counting chamber, dark lined, Bürker	1	631-0920
Cell counting chamber, dark lined, Bürker-Türk	1	630-1543
Cell counting chamber, dark lined, with clamps, Bürker	1	630-1541
Cell counting chamber, dark lined, with clamps, Bürker-Türk	1	631-0922



## Cell counters Marienfeld




### For *in vitro* diagnostic (IVD) applications

These counting chambers are suitable for counting leucocytes and erythrocytes. They are supplied with two cover glasses - Thickness: 0,4 mm.

- Double net ruling: Two counting grids with one division
- Approved for *in vitro* diagnostic (IVD) applications
- Meet the demands of the IVD directive 98/79/EC

Description	Depth (mm)	Pk	Cat. No.
Cell counting chamber, bright lined, without clamps, Neubauer Improved	0,1	1	630-1509
Cell counting chamber, dark lined, without clamps, Fuchs-Rosenthal	0,2	1	630-1510
Cell counting chamber, dark lined, without clamps, Nageotte	0,5	1	630-1511

## Cell counters, Fuchs-Rosenthal



### For counting cells in cerebro-spinal fluid

This ruling is used for counting leukocytes, eosinophils and cells in cerebro-spinal fluid.

- Subdivision: 16 square groups of 1 mm<sup>2</sup> each, consisting of 16 mini squares (0,25×0,25 mm) each

Total area of ruling: 16 mm<sup>2</sup>

Description	Pk	Cat. No.
Cell counter, Fuchs-Rosenthal, depth 0,1 mm	1	631-0930
Cell counter, double, Fuchs-Rosenthal, depth 0,2 mm	1	631-0929
Cell counter, double, with clamps, Fuchs-Rosenthal, depth 0,2 mm	1	631-0924