

Description	Pk	Cat. No.
USC900TH	1	142-6012
USC1200TH	1	142-6010
USC1700TH	1	142-6014
USC2100TH	1	142-6016
USC2600TH	1	142-6018

### Digital models with timer and heater, frequency 45 kHz

Type	Capacity (l)	Max. output power (W)	Tank dimensions WxDxH (mm)	WxDxH (mm)
USC300D	2.75	160	240x137x100	265x162x235
USC500D	4.60	200	240x137x150	265x162x295
USC600D	5.40	260	300x151x150	325x176x295
USC900D	8.60	400	495x130x150	515x150x270
USC1200D	12	400	300x240x200	325x265x335
USC1700D	17	480	327x300x200	352x325x335
USC2100D	21	480	505x300x150	530x325x365
USC2600D	26	660	505x300x200	530x325x365

Description	Pk	Cat. No.
USC300D	1	142-6003
USC500D	1	142-6006
USC600D	1	142-6009
USC900D	1	142-6013
USC1200D	1	142-6011
USC1700D	1	142-6015
USC2100D	1	142-6017
USC2600D	1	142-6019

### Digital models with timer and heater, frequency 132 kHz, for sensitive cleaning of delicate instruments


Type	Capacity (l)	Max. output power (W)	Tank dimensions WxDxH (mm)	WxDxH (mm)
USC300D	2.75	160	240x137x100	265x162x235
USC600D	5.40	260	300x151x150	325x176x295

Description	Pk	Cat. No.
USC300D	1	142-6051
USC600D	1	142-6052

## Accessories for ultrasonic baths



Set comprising of lid and mesh basket



Positioning lids for beakers

Stainless steel and rust-proof. Beakers and test tubes to be cleaned should never lie directly on the bottom of the tank as they will not benefit fully from the cleaning action. Immersion baskets position articles being cleaned correctly in the bath and prevent the bottom of the tank from becoming scratched. Beakers for cleaning smaller parts, or for working with aggressive solutions, can be fitted in positioning covers.

Description	For	Pk	Cat. No.
<b>Set comprising of lid and mesh basket</b>			
	USC100 models	1	142-6048
	USC200 models	1	142-6049
	USC300 models	1	142-6020
	USC500 models	1	142-6021
	USC600 models	1	142-6022
	USC900 models	1	142-6024
	USC1200 models	1	142-6023
	USC1700 models	1	142-6025
	USC2100 models	1	142-6026
	USC2600 models	1	142-6027
<b>Positioning lids for beakers</b>			
Positioning lid for 1 beakers 600 ml	USC200 models	1	142-1704
Coverlid for 2 beakers 250 ml	USC100 models	1	142-6050