



Spencer wells clamps



Stainless steel, rust-proof

For stopping flow in tube connections. Straight with box joint.

- Very good thermal resistance
- Easy cleaning
- Safe, simple operation

Length (mm)	Pk	Cat. No.
125	1	232-0040
152	1	232-0041
178	1	232-0042



Non-Return valves Kartell



PP, translucent

Prevent pressure loss and backflow. Maximum working pressure 2 bar and minimum working pressure 0,07 bar at room temperature (+20 °C).

- Very good chemical resistance
- Very good thermal resistance, autoclavable
- Tapered ends to ensure that tubes are held firmly

For tube I-Ø (mm)	Pk	Cat. No.
7,6 - 10	10	229-3330
9,2 - 12,5	10	229-3331

Non-return valves, with tapered inlet/outlet Bürkle



PP, translucent

Retaining valves with embossed directional arrow or printed in black. With tapered ends. Various uses.

- Resistant to compression max. 4 bar
- Operating temperature up to +40 °C
- With PP ball valve

Type	For tube I-Ø (mm)	Pk	Cat. No.
With embossed directional arrow	6 - 8	1	229-0134
With embossed directional arrow	8 - 10	1	229-0135
With embossed directional arrow	10 - 14	1	229-0136
With black printed directional arrow	8 - 10	1	229-0137

ALL YOU NEED FOR FOOD ANALYSIS

Specialist publication for all those working in food science

Request your copy from your local VWR sales office or vwr.com