

Continued from previous page

Application	Manual Cleaning					Machine cleaning							
Phosphate analysis equip.			x										
Pipettes		x		x		x		x					
Plaster residues	x	x	x		x								
Plastic equipment						x		x				x	
Precision equipment		x		x									
Protein residues	x		x	x			x						
Proteins	x		x	x					x				
Quartz equipment		x											
Rubber		x											
Saliva	x		x	x									
Silicones	x		x						x		x		
Spectacle lenses		x											
Stainless steel	x	x	x	x	x				x				
Thin film plates	x		x		x	x							
Tiled surfaces	x		x		x				x				
Tough residues				x			x						
Ultrasound	x	x	x										
Wax	x		x	x									
Zinc		x											

Description	pH value	Temp. range (°C)	Pk	Cat. No.
<b>Manual cleaning - immersion cleaning</b>				
LABWASH® Premium Classic, mildly alkaline concentrate for manual cleaning of laboratory glassware	12,2	20...85	2,5 l	84545.320
	12,2	20...85	5 l	84545.360
	12,2	20...85	10 l	84545.410
LABWASH® Premium Extra, concentrate for the manual cleaning of laboratory glassware and precision components	11,8	10...40	2,5 l	84546.320
	11,8	10...40	5 l	84546.360
	11,8	10...40	10 l	84546.410
LABWASH® Premium Extra PF, phosphate-free concentrate for manual cleaning of laboratory glassware	12,7	10...70	2,5 l	84554.320
	12,7	10...70	5 l	84554.360
	12,7	10...70	10 l	84554.410
LABWASH® Premium Neutral, pH-neutral concentrate for the cleaning of laboratory glassware and precision components	7,4	10...60	2,5 l	84547.320
	7,4	10...60	5 l	84547.360
	7,4	10...60	10 l	84547.410
<b>Automated cleaning - laboratory washing machines</b>				
LABWASH® Premium Powder, mildly alkaline cleaning powder for laboratory washing machines	12,1	10...60	1 kg	84548.290
	12,1	10...60	10 kg	84548.410
LABWASH® Premium Powder PF, phosphate-free, alkaline powder cleaner for laboratory washing machines	13	0...95	1 kg	84557.290
	13	0...95	10 kg	84557.410
LABWASH® Premium Pure, concentrate for soak and spray cleaning of laboratory glassware	11,8	0...95	2,5 l	84555.320
	11,8	0...95	5 l	84555.360
	11,8	0...95	10 l	84555.410
LABWASH® PREMIUM Universel, highly alkaline liquid cleaning concentrate for laboratory washing machines	13,3	10...95	10 l	84717.410
	13,3	10...95	20 l	84717.415
LABWASH® Premium Alkamatic, mildly alkaline, liquid cleaning concentrate for washing machines	12	20...95	2,5 l	84549.320
	12	20...95	5 l	84549.360
	12	20...95	10 l	84549.410
LABWASH® Premium Alkamatic LA, highly alkaline, intensive cleaner for laboratory washing machines	12,5	20...95	2,5 l	84556.320
	12,5	20...95	5 l	84556.360
	12,5	20...95	10 l	84556.410
LABWASH® Premium Alkamatic PF, alkaline, phosphate-free, chlorine-free, liquid cleaning concentrate for washing machines for heavy contamination	11,6	10...95	2,5 l	84550.320
	11,6	10...95	5 l	84550.360
	11,6	10...95	10 l	84550.410
LABWASH® PREMIUM Neutraliz, acetic acid based neutraliser and pre-cleaner for laboratory washing machines	11,6	10...95	20 l	84550.440
	2,7	40...85	10 l	84718.410
	2,7	40...85	20 l	84718.415
LABWASH® Premium Acidrinse P, neutraliser and pre-cleaner for washing machines	2,1	40...85	2,5 l	84552.320
	2,1	40...85	5 l	84552.360
	2,1	40...85	10 l	84552.410
LABWASH® Premium Acidrinse C, phosphate-free neutraliser for washing machines	2,7	40...85	2,5 l	84551.320
	2,7	40...85	5 l	84551.360
	2,7	40...85	10 l	84551.410
LABWASH® Premium Acidrinse C, phosphate-free neutraliser for washing machines	2,7	40...85	20 l	84551.440
	2,2	0...60	2,5 l	84553.320
	2,2	0...60	5 l	84553.360
LABWASH® Premium Easyrinse, acidic, ecological rinsing concentrate with shiny effect	2,2	0...60	10 l	84553.410

pH value at 1% concentration