

## Staining jars Kartell



### PMP (TPX), transparent

Staining jar similar to Hellendahl design for 16 standard microscope slides 76×26 mm (in pairs, back to back).

- Extension and cover

Description	WxDxH (mm)	Pk	Cat. No.
Staining jar, Hellendahl type	57×57×90	1	631-9130

## Staining jars according to Hellendahl Wheaton



### Soda-lime glass

These staining jars hold up to eight single standard slides (76×26 mm) positioned vertically or 16 slides back to back.

- Wide top for convenient transfer
- For staining slides that are inscribed on one end

Description	WxDxH (mm)	Pk	Cat. No.
Staining jar acc. to Hellendahl	76×25×75	6	720-0711

## Staining jars according to Hellendahl Wheaton



### Soda-lime glass

Staining jar holds eight individual 76×26 mm standard slides or 16 slides back to back.

- Easy to use due to wide opening
- For staining or use as a TLC developing chamber

Description	WxDxH (mm)	Pk	Cat. No.
Staining jar acc. to Hellendahl	42×26×85	6	720-0710

## Staining jars



### Soda-lime glass

For 10 standard microscopy slides (76×26 mm).

Description	Pk	Cat. No.
Staining jar according to Coplin	1	631-9331