

## Chromatography columns with Rotaflo® stopcock



PYREX®, borosilicate glass

- With PTFE Rotaflo® stopcock for fine control of eluent flow
- Chemically resistant glass sintered disc - porosity 0
- Ground glass socket accepts Quickfit® reservoirs

Socket size (NS)	I-ØxL (mm)	Pk	Cat. No.
14/23	10x100	1	201-0722
14/23	10x200	1	201-0723
14/23	10x300	1	201-0724
14/23	10x400	1	201-0725
19/26	20x200	1	201-0727
19/26	20x300	1	201-0728
19/26	20x400	1	201-0729
19/26	20x500	1	201-0730
24/29	26x400	1	201-0731
24/29	26x500	1	201-0732
29/32	30x400	1	201-0733
29/32	30x500	1	201-0734
34/35	30x400	1	201-0735
34/35	30x500	1	201-0736
40/38	40x400	1	201-0737
40/38	40x500	1	201-0738



## Chromatography columns with Rotaflo® stopcock and ground cone



PYREX®, borosilicate glass

- With glass sintered disc (porosity 0) to support column packing
- GP Rotaflo® stopcock with PTFE key for fine control of eluent flow
- Stem within lower cone allows mounting on a Büchner flask for vacuum flash chromatography

Cone size (NS)	Socket size (NS)	Total length (mm)	I-ØxL (mm)	Pk	Cat. No.
19/26	14/23	270	10x100	1	201-8297
19/26	14/23	470	10x300	1	201-0726



VWR.com in your pocket

### VWR2GO

VWR2Go provides the essential features from our web-based shop system on mobile phones and tablets

Real time, intuitive searching and ordering of over one million items

Download today!