

Continued from previous page

- Antimicrobial surface - silver ions (Ag) create a surface which inhibits growth of bacteria, microbes and fungus
- Advanced Volume Gearing mechanism (ACG) - self supportive modular mechanism increasing accuracy, precision and durability, with light force and smaller increments for volume adjustment
- Set-and-forget pipetting button - easy set and forget lock mechanism to adjust the volume, pipetting button rotates freely when not in locking position
- Superior handle ergonomics, including adjustable finger rest (120°) for optimal pipetting position, very low force pipetting movement
- Super blow-out piston for volumes of 50 µl and below, which ensures the delivery of micro-size drops

Ordering information: Supplied with one opening and calibration key, one maintenance set (incl. grease and O-rings), one calibration report, three-year warranty certificate (extended to five years with web registration) and user's guide (English, French, German, Spanish).

GLP-Kits supplied with three Finnpipe[®] F1 (GLP-Kit 1 and 3) or four Finnpipe[®] F1 (GLP-Kit 2), one stand for 6 pipettes, three racks of Finn[®]tip Flex tips and 15 disposable reagent reservoirs.

Capacity (µl)	Division (µl)	Accuracy (%)	Imprecision (%)	Colour code	For Finn [®] tip tip(s)	Pk	Cat. No.
Fixed volume							
1*	-	±4,00	<4,00	—	Flex 10, 10, 20/50 Micro	1	613-2742
5*	-	±1,40	<1,40	—	Flex 10, 10, 20/50 Micro	1	613-2743
10	-	±0,90	<0,80	—	Flex 200, 250 Univ., 200 Ext.	1	613-2744
20	-	±0,70	<0,50	—	Flex 200, 250 Univ., 200 Ext.	1	613-1478
25	-	±0,60	<0,50	—	Flex 200, 250 Univ., 200 Ext.	1	613-2745
50	-	±0,60	<0,40	—	Flex 200, 250 Univ., 200 Ext.	1	613-2746
100	-	±0,40	<0,30	—	Flex 200, 250 Univ., 200 Ext.	1	613-2747
200	-	±0,40	<0,30	—	Flex 200, 250 Univ., 200 Ext.	1	613-1479
250	-	±0,40	<0,30	—	Flex 1000, 1000, 1000 Ext.	1	613-2748
500	-	±0,30	<0,30	—	Flex 1000, 1000, 1000 Ext.	1	613-2749
1000	-	±0,30	<0,30	—	Flex 1000, 1000, 1000 Ext.	1	613-2750
2000	-	±0,30	<0,20	—	5 ml	1	613-2751
3000	-	±0,30	<0,20	—	5 ml	1	613-2752
5000	-	±0,30	<0,20	—	5 ml	1	613-2753
10000	-	±0,30	<0,20	—	10 ml, Flex 10 ml Ext.	1	613-2754
Variable volume							
0,2 - 2*	0,002	±2,50 - ±12,00	<2,00 - <0,00	Pink	Flex 10, 10, 20/50 Micro	1	613-2730
0,5 - 5*	0,01	±1,50 - ±6,00	<1,00 - <5,00	Pink	Flex 10, 10, 20/50 Micro	1	613-2731
1 - 10*	0,02	±1,00 - ±2,50	<0,50 - <2,00	Pink	Flex 10, 10, 20/50 Micro	1	613-2732
1 - 10	0,02	±1,00 - ±3,50	<0,80 - <3,00	Yellow	Flex 200, 250 Univ.	1	613-2733
2 - 20*	0,02	±1,00 - ±3,00	<0,40 - <2,50	Turquoise	20/50 Micro	1	613-2734
2 - 20	0,02	±1,00 - ±3,00	<0,40 - <2,50	Yellow	Flex 200, 250 Univ.	1	613-2735
5 - 50*	0,1	±0,60 - ±3,00	<0,30 - <2,50	Turquoise	50 Micro	1	613-1484
5 - 50	0,1	±0,60 - ±3,00	<0,30 - <2,50	Yellow	Flex 200, 250 Univ.	1	613-1485
10 - 100	0,2	±0,80 - ±3,00	<0,20 - <1,00	Yellow	Flex 200/300, 250 Univ., 200 Ext.	1	613-2736
20 - 200	0,2	±0,60 - ±1,80	<0,20 - <0,70	Yellow	Flex 200/300, 250 Univ., 200 Ext.	1	613-2737
30 - 300	1	±0,60 - ±1,60	<0,20 - <0,60	Orange	Flex 300, 300	1	613-2738
100 - 1000	1	±0,60 - ±1,00	<0,20 - <0,60	Blue	Flex 1000, 1000, 1000 Ext.	1	613-2739
500 - 5000	10	±0,50 - ±2,00	<0,20 - <0,80	Green	5 ml	1	613-2740
1000 - 10000	20	±0,50 - ±2,00	<0,20 - <0,80	Red	10 ml, Flex 10 ml Ext.	1	613-2741
Description		For			Pk		Cat. No.
Pipette stands for Finnpipe[®]							
F-Stand for 6 F1/F2/F3 Finnpipe [®]		F1/F2/F3 Finnpipe [®]			1		613-2765

* Micro