



Laboratory stands



Tripod stands (A): Cast iron tripod stand, painted with hammer finish, steel support rod, nickel plated

- Threaded hole
- Support rod (three lengths available with different diameters)
- Robust, corrosion resistant

Support stands (B): Cast iron stand base, painted with hammer finish, steel support rod

With eccentrically attached support rod.

- Firm, non-slip stand
- Available in three different sizes
- Robust, corrosion resistant

Support stands (C): Cast iron stand base, painted with hammer finish, steel support rod, nickel plated


With centrally attached support rod.

- Firm, non-slip stand
- Available in three different sizes
- Robust, corrosion resistant

Description	Length (mm)	Width (mm)	Arm length (mm)	Pk	Cat. No.
Tripod stands (A)					
Tripod stand (A) with rod (500 mm length, 8 mm Ø)	-	-	115	1	241-0090
Tripod stand (A) with rod (800 mm length, 12 mm Ø)	-	-	145	1	241-0091
Tripod stand (A) with rod (1000 mm length, 16 mm Ø)	-	-	175	1	241-0092
Support stands (B)					
Support stand (B) with rod (500 mm length, 8 mm Ø)	225	140	-	1	241-0093
Support stand (B) with rod (800 mm length, 12 mm Ø)	250	160	-	1	241-0094
Support stand (B) with rod (1000 mm length, 16 mm Ø)	270	180	-	1	241-0095
Support stands (C)					
Support stand (C) with rod (500 mm length, 8 mm Ø)	300	140	-	1	241-0099
Support stand (C) with rod (800 mm length, 12 mm Ø)	340	160	-	1	241-0100
Support stand (C) with rod (1000 mm length, 16 mm Ø)	380	180	-	1	241-0101



Stand with holder for separating funnels



PU coated cast iron stand 200×450 mm, steel rod 500×12 mm, bosshead and holder for two funnels.

- Height of holder can be adjusted
- Holder measures 300×100 mm
- Apertures for funnels measure 65 mm - centre distance 200 mm

Description	L×W (mm)	Pk	Cat. No.
Stand with funnel holder for two funnels	200×450	1	241-0046



Stand with holder for funnels



PU coated cast iron stand (200×450 mm), steel rod (500×12 mm), bosshead and holder for four funnels.

- Height of holder can be adjusted
- Holder measures 400×80 mm
- Apertures for funnels measure 100 mm in diameter and recesses in front of the apertures for funnels have a width of 40 mm

Description	L×W (mm)	Pk	Cat. No.
Stand with funnel holder for four funnels	200×450	1	241-0045