

Centrifuges 5702, 5702 R and 5702 RH

Eppendorf



The 5702 range of centrifuges is ideal for medical routine, clinical and biotechnology research. They can be used with 1.1 ml tubes for blood samples and with Falcon® tubes. Rotors, buckets and lids are autoclavable. They have induction motors and a timer which can be set for up to 90 minutes or for continuous operation.



- User friendly control panel with a small number of clearly labelled function keys/knobs and easy to read LCD display; speed can be set from 100 rpm to a maximum of 4400 rpm
- Automatic motorised locking of the lid
- Slow startup and braking through the SOFT function
- Stainless steel rotor chamber (rust proof)
- 5702 RH can consistently maintain special temperatures above room temperature, temperature-sensitive cells or cell cultures can be centrifuged at defined temperatures, particularly at 37 °C

Model	5702	5702 R	5702 RH
Max. Speed (t/min)		4400	
Max. RCF (g)		3000	
Capacity		4x90 ml	
Temperature (°C)		-9...+40	-9...+42
Acceleration/ braking times (s)		< 25/25	
Sound level (dB)	< 58	< 54	
Power output (W)	200	380	
WxDxH (mm)	320x400x240	380x580x260	
Weight (kg)	20	36	

Description	Pk	Cat. No.
Centrifuge 5702, without rotor, UK-plug	1	521-0125
Centrifuge 5702 R, refrigerated, without rotor, UK-plug	1	521-0134
Centrifuge 5702 RH refrigerated and heatable, without rotor	1	521-0593

Accessories

Description	Pk	Cat. No.
Adapter 1 x 15 ml conical (Ø 17.2 mm), 3000 g	2	521-0037
Adapter 1 x 25 ml (Ø 25 mm), 2900 g	2	521-0029
Adapter 1 x 50 ml conical (Ø 30 mm), 2900 g	2	521-0038
Adapter 1 x 85 ml (Ø 38 mm), 2900 g	2	521-0030
Adapter 2x15 ml conical (Ø 17.2 mm), 3000 g (use without cap)	2	521-0035
Adapter 4 x 2.6/7 ml (Ø 13.5 mm), 2850 g	2	521-0031
Adapter 4 x 4/10 ml (Ø 16 mm), 2850 g	2	521-0039
Adapter 4 x 9/15 ml (Ø 17.5 mm), 2800 g	2	521-0036
Adapter 5 x 1.1/1.4 ml (Ø 8.5 mm), 2850 g	2	521-0040
Adapter 5 x 2/7 ml (Ø 12.5 mm), 2850 g	2	521-0041
Adapter 6 x 15 ml (Ø 17.5 mm), 2750 g	2	522-1153
Aerosol-tight lid for 85 ml round bucket	2	521-0033
Fixed-angle rotor F-35-30-17 10 x 15 ml, with 10 Falcon® adapters and 30 rubber mats, 2750 g	1	521-0025
Fixed-angle rotor F-35-30-17 30 x 15 ml, with 20 Falcon® adapters and 30 rubber mats, 2750 g	1	521-0024
Fixed-angle rotor F-45-18-17-Cryo 18 x 1.0/2.0 ml, without lid, 1840 g	1	521-0597
Fixed-angle rotor F-45-24-11 24 x 1.5/2.0 ml, without lid, 1770 g	1	521-0595
Rectangular buckets 90 ml	4	522-1132
Swing-bucket rotor A 4-36 without round buckets, 3000 g	1	521-0133
Swing-bucket rotor A 4-38 with 4 round buckets 85 ml, 3000 g (adapters not included)	1	521-0032
Swing-bucket rotor A-8-17, 8x15 ml with 8 inserts, 2800 g	1	521-0021