



Eyewash solution Plum



500 ml



200 ml

Sterile sodium chloride solution (0,9%) in an eyewash bottle with eye cup and dust cap. Helps prevent foreign bodies (such as metal and wood chips) from sticking to the eye and causing serious eye injury.

- Ergonomically formed eye cup fits the eye and enables gentle and even rinsing
- 200 ml bottle particularly suitable for mobile use, such as in tool boxes or first aid kits
- Can be used to refill Plum eyewash stations
- No maintenance and cleaning of the bottles required

Rinsing time: Approx. 2 minutes per 200 ml, 5 minutes per 500 ml

Shelf life: 3 years

According to DIN EN 15154-4

Description	Pk	Cat. No.
Eyewash solution, 200 ml	1	118-1703
Eyewash solution, 500 ml	1	118-1700

Description	Pk	Cat. No.
Accessories		
Belt bag for 200 ml eyewash bottles	1	118-1704
Green wall-mountable holder for 200 ml and 500 ml eyewash bottles	1	118-1711

Eyewash cabinets and stations



Two types of emergency eyewash cabinets supplied with 2x500 ml bottles of sterile saline solution and a choice of two or four eye pads. The wall-mounted eyewash station is completed with 2x500 ml sterile saline solution, displayed instructions and a fitted mirror.

- Dust and dirt proof
- Choice of two cabinets and one emergency eyewash station
- Ideal for use in a workplace with fewer users

Ordering information: Spare solution bottles should be ordered separately.

Description	Pk	Cat. No.
Eyewash cabinet with 4 sterile eye pads and 2 x 500 ml sterile eyewash	1 KIT	118-0115
Eyewash cabinet with 2 sterile eye pads and 2 x 500 ml sterile eyewash	1 KIT	118-0117
Eyewash station complete, with sterile saline solution bottles, mirror and instructions	1 KIT	118-0118

Description	Pk	Cat. No.
Spare parts		
Eye wash solution, refill 500 ml	1	118-0114

Dispenser for eyewash pods



Wall-mounted eyewash pod dispenser with integrated mirror for easier self-administration.

- Clear instructions for use printed on front of cabinet
- Dust and dirt proof
- Will hold up to 30 pods (supplied empty)

Ordering information: Sterile pods (118-0113) need to be ordered separately.

Description	Pk	Cat. No.
Dispenser without pods	1	118-0116

Description	Pk	Cat. No.
Sterile pods		
Eye wash pods 20 ml	25	118-0113