

Continued from previous page

- User friendly control panel with a small number of clearly labelled function keys and easy to read LCD display; speed can be set from 100 min⁻¹ to a maximum of 4400 min⁻¹
- Automatic, motorised lid locking
- Slow startup and braking through the SOFT function
- Stainless steel rotor chamber (rust-proof)
- Heated and refrigerated centrifuge 5702 RH allows centrifugation under biological conditions to improve cell viability for more accurate results in downstream applications (e.g. cell cultivation)

Model	5702	5702 R	5702 RH
Max. capacity		4×90 ml	
Max. speed (min ⁻¹)		4400	
Max. RCF (×g)		3000	
Temperature (°C)	-	-9...+40	-9...+42
Acceleration/braking time (s)		<25/25	
Noise level (dB(A))	<51		<45
Power consumption (W)	200		380
W×D×H (mm)	320×400×240		380×580×260
Weight (kg)	20		36

Description	Pk	Cat. No.
Centrifuge 5702, without rotor, UK-plug	1	521-0125
Centrifuge 5702 R, refrigerated, without rotor, UK-plug	1	521-0134
Centrifuge 5702 RH, refrigerated and heated, without rotor, UK-plug	1	521-0594

Description	Pk	Cat. No.
Swing-out rotors		
Swing out rotor A-8-17, 8×15 ml with 8 inserts, 2800 ×g	1	521-0021
Swing out rotor A-4-38 with 4×85 ml round buckets, 3000 ×g (adapters sold separately)	1	521-0032
Swing out rotor A-4-38 without round buckets, 3000 ×g	1	521-0133
Fixed-angle rotors		
Fixed angle rotor F-35-30-17, 30×15 ml, with 30 conical tube adapters and 30 rubber mats, 2750 ×g	1	521-0024
Fixed angle rotor F-35-30-17, 10×15 ml, with 10 conical tube adapters and 10 rubber mats, 2750 ×g	1	521-0025
Fixed angle rotor F-45-24-11, 24×1,5/2,0 ml, without lid, 1770 ×g	1	521-0595
Fixed angle rotor F-45-18-17-Cryo, 18×1,0/2,0 ml, without lid, 1840 ×g	1	521-0597
Adapter		
Rectangular buckets, 90 ml (adapters sold separately) for A-4-38 rotor	4	522-1132
Aerosol-tight lids for 85 ml round buckets	2	521-0033

IND

Centrifuge, bench top, ROTOFIX 32 A Hettich



A rugged, versatile and indispensable centrifuge for routine laboratory tasks in doctors' practices and small hospital laboratories, especially for tasks in clinical chemistry and cytology as well as for the preparation of samples in industrial and research laboratories.

- Metal housing and lid, stainless steel chamber, practical lid locking device, viewing port in the lid
- Easy to use foil keypad for fast, precise and easy entry of the parameters via selector keys
- Pulse key for short centrifuging, RCF key enabling indication of the RCF during centrifugation
- Safety features: Lid locking and holding, imbalance switch off, automatic rotor recognition, emergency lid lock release

Ordering information: Rotors and adapters available separately.

Model	ROTOFIX 32 A
Max. capacity - swing out rotor	4×100 ml / 32×15 ml
Max. speed - fixed angle rotor (min ⁻¹)	6000
Max. RCF - fixed angle rotor (×g)	4226
Timer	1 - 99 min, continuous run, short cycle mode (pulse key)
Power supply	208 - 240 V / 50 - 60 Hz
W×D×H (mm)	366×430×257
Weight (kg)	23

Description	Pk	Cat. No.
Bench top centrifuge ROTOFIX 32 A, without rotor	1	521-0902

Description	Pk	Cat. No.
Swing out rotors		
Swing out rotor 1324, 4×100 ml, 4500 min ⁻¹ , 3328 ×g, without buckets	1	521-1303
Swing out rotor 1624, 4×100 ml, 4000 min ⁻¹ , 2665 ×g, to be completed with buckets and adapters	1	523-1011

Continued on next page