

Gallium standard solution, 10,000 mg/l Ga in dil. nitric acid (from Ga) ARISTAR® standard for ICP

Ga in HNO₃

Traceable to SRM from NIST, tested in ISO Guide 34 / ISO17025 accredited laboratory.

Cat. No.	Pk	Pack type
455392G	100 ml	Plastic bottle

Supplied with certificate of analysis.

Gallium standard solution, 1,000 mg/l Ga in dil. nitric acid

Warning

H319 H315

P280 P302+P352 P305+P351+P338

CAS 7440-55-3

EINECS: 231-163-8

UN: 3264

ADR 8,III

Ga

M.W. 69.723 g/mol

Density: 1.014 g/cm³ (20 °C)

Storage Temperature: Ambient temperature



Gallium standard solution, 1,000 mg/l Ga in dil. nitric acid (from Ga) ARISTAR® standard for ICP

Ga in HNO₃

Traceable to SRM from NIST, tested in ISO Guide 34 / ISO17025 accredited laboratory.

Cat. No.	Pk	Pack type
455382E	100 ml	Plastic bottle
455384G	500 ml	Plastic bottle

Supplied with certificate of analysis.

NEW Gallium standard solution, 1,000 mg/l Ga in dil. nitric acid AVS TITRINORM® standard for AAS

Cat. No.	Pk	Pack type
86678.180	100 ml	Plastic bottle
86678.260	500 ml	Plastic bottle

Supplied with certificate of analysis.

Gardner Colour Standard (2-16)

Storage Temperature: Ambient temperature

NEW Gardner Colour Standard 2

Cat. No.	Pk	Pack type
84815.260	500 ml	Glass bottle

NEW Gardner Colour Standard 4

Cat. No.	Pk	Pack type
84816.260	500 ml	Glass bottle

NEW Gardner Colour Standard 6

Cat. No.	Pk	Pack type
84817.260	500 ml	Glass bottle

NEW Gardner Colour Standard 8

Cat. No.	Pk	Pack type
84818.260	500 ml	Glass bottle

NEW Gardner Colour Standard 10

Cat. No.	Pk	Pack type
84819.260	500 ml	Glass bottle

NEW Gardner Colour Standard 12

Cat. No.	Pk	Pack type
84820.260	500 ml	Glass bottle

NEW Gardner Colour Standard 14

Cat. No.	Pk	Pack type
84821.260	500 ml	Glass bottle

NEW Gardner Colour Standard 16

Cat. No.	Pk	Pack type
84822.260	500 ml	Glass bottle

E-BUSINESS

More than just a webshop for all of your laboratory needs

Please visit us at vwr.com

AVS® TITRINORM® READY TO USE SOLUTIONS

- Volumetric solutions
- pH buffers
- AAS standard solutionsAv