Swab rinse kit



Swab rinse kits (SRK) for microbiological examination of surfaces. The isotonic salt solution contains substances to neutralise and inactivate disinfectants and sanitising agents. The rinse solution is used to moisten the sterile swab before sampling and to transport the swab after collection.



· Can be used for microbial detection or quantitative analysis

Short SRK swabs (rayon tip, attached to the tube cap)

Each box contains 50 screw cap tubes. Swab shafts are attached to the tube caps.

Description	Packed	Length (mm)	Pk	Cat. No.	
SRK rinse solution, 10 ml	single wrapped in plastic peel pouch	80	50	710-0477	
SRK rinse solution, 10 ml	50 tubes/ box	80	50	710-0478	

Long SRK swabs (14.6 cm - separate from cap)

Each box contains 25 pouches (each containing a tube and swab). Swabs are separate from cap.

Description	Packed	Length (mm)	Pk	Cat. No.	
SRK rinse solution, rayon, 10 ml	tube and applicator in plastic peel pouch	146	25	710-0472	
SRK rinse solution, alginate, 10 ml	tube and applicator in plastic peel pouch	146	25	710-0473	
SRK rinse solution, rayon, 2.5 ml	tube and applicator in plastic peel pouch	146	25	710-0431	
SRK rinse solution, alginate, 2.5 ml	tube and applicator in plastic peel pouch	146	25	710-0476	

Large SRK swabs (19.50 cm long - separate from cap)

Each box contains 25 pouches (each containing a tube and swab). Swabs are separate from cap.

Description	Packed	Length (mm)	Pk	Cat. No.	
SRK rinse solution, fibre tip rayon, 10 ml	tube and applicator in plastic peel pouch	195	25	710-0474	
SSK square sampling kit, fibre tip rayon, 10 ml	tube, large applicator & sterile 10×10 cm template mask in sterile peel pouch	195	1 Kit	710-0475	

SRK foam spatulas

Small spatula attached to cap. Large foam spatula, 22 cm long (612-2678).

Description	Packed	Pk	Cat. No.	
SRK rinse solution, 10 ml	in plastic peel pouch	50	612-2679	
SRK without medium, short	in plastic peel pouch	50	612-2680	
SRK, without medium, large	in plastic peel pouch	120	612-2678	

SRK 471

Kit for bacteriological sampling of carcasses. Designed in accordance with EC directive 471/2001/EC. Eight plain swabs in tubes with blue shafts. Packed together with one tube of SRK rinse solution in a zip seal bag.

Description	Packed	Lenath (mm)	Pk	Cat. No.	
SRK rinse solution, ravon, 10 ml	6×25 packs per case	145	150	710-0479	

Swabs, Steriswabs®



Blue plastic shafted swabs, with 45 mm breakpoint, dry or pre-moistened with sterile diluent, designed for effective monitoring of sterility and contamination in clinically clean or sterile work areas.



www vwr com

- Available in labelled tube or peel pouches
- · Triple wrapped and radiation sterilised

Description	Pk	Cat. No.	
Swab, blue plastic shaft, labelled tube	250	720-0142	