



Absorbent pad kits



One-handed dispensing of cellulose absorbent pads Ø 45.5 mm.

Absorbent pads are ideal for absorbing broth media to culture colonies in accordance with the membrane filter technique.

- Enables user to dispense a clean cellulose pad into a Petri dish whenever needed without touching the pad
- Dispenser kit holds one tube of one hundred absorbent pads (ten tubes included)
- Available non sterile or gamma irradiated (no ethylene oxide residuals to impede microbial growth)

Description	Pk	Cat. No.
Absorbent pad kit, sterile	1.000	516-9012
Absorbent pad kit, non sterile	1.000	516-7850

Petri dish carrying rack, Poxygrid® Bel-Art Products



Epoxy coated steel wire rack

- Easily loaded and emptied from the top
- For plastic Petri dishes Ø 90 and 100 mm (not suitable for glass Petri dishes)

WxDxH: 108x108x356 mm (1 stack, for 20 dishes)

WxDxH: 343x105x305 mm (3 stacks, for 42 dishes)

Description	Pk	Cat. No.
Petri dish carrying rack, 1 stack	1	391-2011
Petri dish carrying rack, 3 stacks	1	391-2013

VWR.COM TALKING YOUR LANGUAGE

The local website with global reach