

Water Purification

Water Demineralisation

Deionisers, Vision 125 and 250

ELGA LabWater



Vision deionisers provide quality water as and when it is required. The need for elaborate and costly storage is eliminated. Just turn on the mains water tap and take off the water instantly either through the top dispensing spout or the optional hand-held dispensing gun on a flexible lead. The resin changes colour when the cartridge is exhausted.

- Windows enable users to monitor resin colour and cartridge status
- High quality water instantly available
- Simple to install, wall mounting

| Model | Vision 125 | Vision 250 |
|-----------------------------------|------------|------------|
| Max. flow rate (l/h) | 60 | 60 |
| Max. operating pressure (psi/bar) | 90/6 | 90/6 |
| Typical water quality (µS/cm) | 0.1 to 10 | 0.1 to 10 |
| WxDxH (mm) | 92x110x385 | 92x110x660 |
| Weight (kg) | 2.5 | 4.2 |

Delivery Information: Supplied with dispensing gun.

| Description | Pk | Cat. No. |
|--------------------------------------|----|----------|
| Deioniser, Vision 125, wall-mounted | 1 | 171-0538 |
| Deioniser, Vision 250, wall-mounted | 1 | 171-0539 |
| Replacement cartridge for Vision 125 | 1 | 171-0536 |
| Replacement cartridge for Vision 250 | 1 | 171-0537 |
| Gun Dispenser for Vision deionisers | 1 | 171-0520 |

Deionisers, Micromeg

ELGA LabWater



This wall mounted disposable deioniser produces water up to a quality of 10M.-cm at flows of up to 90 l/hr. Designed for users with a requirement for low volumes deionised water, the Micromeg requires no electrical input as its water quality meter is battery operated.

- Disposable cartridge
- New cartridge can be installed in seconds

| | |
|-----------------------------------|-------------|
| Max. flow (l/h) | 90 |
| Max. operating pressure (psi/bar) | 50/3.5 |
| WxDxH (mm) | 200x230x550 |
| Weight (kg) | 6.5 |

| Description | Pk | Cat. No. |
|--|----|----------|
| Deioniser, Micromeg | 1 | 171-0484 |
| Deioniser, Micromeg with exchangeable cartridge | 1 | 171-0485 |
| Deioniser, Micromeg with disposable cartridge | 4 | 171-0486 |
| Deioniser, Micromeg, strong base with transit case | 1 | 171-0487 |