

Membrane dispenser for microbiological control, Membrane-Butler Whatman (GE Healthcare)



The Membrane-Butler dispenses individual sterile MicroPlus and ME membrane filters (type 'STL'). The sterile packaging is inserted into the roller system and the system is ready to be used.

- Compact dimensions for portable use
- Stainless steel housing, ideal for use on sterile benches
- Cross-contamination risks are minimised

Description	Pk	Cat. No.
Membrane-Butler membrane dispenser, manual	1	511-5500

Membrane dispenser for microbiological filtration, Microsart® e.motion Sartorius



The dispenser releases membrane filters from their packaging. Two modes of fully automatic operation: Touch-button or hands-free. An optical sensor detects approaching tweezers.

- Compact design, easy to clean
- Easy insertion of the filter band
- Low weight

Membranes Microsart® e.motion

Cellulose nitrate, sterile

The membrane filters are sterile-sealed, without protective paper on top of each filter, in a specially designed individual package in a fanfold format.

Description	Pk	Cat. No.
Automatic dispenser for Microsart® e.motion	1	511-0009

Description	Pore size (µm)	Ø (mm)	Pk	Cat. No.
Membranes Microsart® e.motion				
Green with dark green grid	0,45	47	300	514-0078
Grey with white grid	0,45	47	300	514-0077
Grey with white grid	0,65	47	300	514-0058
White with black grid	0,45	47	300	514-0076
White with black grid	0,45 high flow	47	300	514-7015
White with green grid	0,45	47	300	514-0079
White with green grid	0,45 high flow	47	300	514-0057



Membrane filter dispenser, Sentino™



Filter dispensers are ideal for high throughput microbiology testing laboratories looking for more ergonomic, faster ways to handle individual membranes. The Sentino™ filter dispenser is designed to aseptically dispense one disc of membrane at the touch of a button. The Sentino™ filter dispenser can be used for:

- Drinking and wastewater analysis
- Environmental and recreational water monitoring
- Personal care product manufacturing
- Beverage monitoring
- Pharmaceutical and industrial water system monitoring
- Pharmaceutical raw materials, intermediates and final product

Delivery information: Supplied with power transformer, EU-plug, UK-plug and user guide CD.

Description	Pk	Cat. No.
Membrane filter dispenser, Sentino™, EU/UK-plug	1	516-0665

Type	Pore size (µm)	Pk	Cat. No.
Dispenser pack refills for membrane filter dispenser, Sentino™			
Supor®	0,2	1.000	514-0629
GN-6 Metricel®	0,45	1.000	514-0628
Metricel® Black	0,45	1.000	514-0630
Metricel® Black	0,8	1.000	514-0631