




## Filter tips, ART® reload system

### Thermo Scientific ART™

This sterile barrier tip reload system combines individually packaged blister tubs made of recycled plastic and new hinged racks. The push down blister tub easily inserts the ART® barrier tips into the hinged racks providing stackable space-saving racks that offer easy, one-handed use. Empty hinged racks allow for easy conversion to the ART® Reload system.

- Environmentally friendly packaging
- ART® self-sealing barrier preventing cross-contamination
- Tips are certified DNase, RNase and pyrogen-free

|                                                                                          | Type                       | Capacity (µl) | Sterile | Packed                            | Pk        | Cat. No.        |
|------------------------------------------------------------------------------------------|----------------------------|---------------|---------|-----------------------------------|-----------|-----------------|
| <b>ART® 10</b>                                                                           |                            |               |         |                                   |           |                 |
|          | MicroPoint™, Low Retention | 10            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1180        |
|                                                                                          | MicroPoint™                | 10            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1147        |
| <b>ART® 10 REACH</b>                                                                     |                            |               |         |                                   |           |                 |
|          | Extended, Low Retention    | 10            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1178        |
|                                                                                          | Extended                   | 10            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1148        |
| <b>ART® 20E</b>                                                                          |                            |               |         |                                   |           |                 |
|          | Ultra-Micro, Low Retention | 10            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1176        |
|                                                                                          | Ultra-Micro                | 10            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1149        |
| <b>ART® 20P</b>                                                                          |                            |               |         |                                   |           |                 |
|          | Low Retention              | 20            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1174        |
|                                                                                          | Standard                   | 20            | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1150        |
| <b>ART® 100</b>                                                                          |                            |               |         |                                   |           |                 |
|          | Standard                   | 100           | +       | 10 racks, 96 each (Hinged rack)   | 960       | 732-1143        |
|                                                                                          | Standard                   | 100           | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1151        |
| <b>ART® 100E</b>                                                                         |                            |               |         |                                   |           |                 |
|        | MicroPoint™, Low Retention | 100           | +       | 10 racks, 96 each (Hinged rack)   | 960       | 732-1168        |
|                                                                                          | MicroPoint™, Low Retention | 100           | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1167        |
|                                                                                          | MicroPoint™                | 100           | +       | 10 racks, 96 each (Hinged rack)   | 960       | 732-1144        |
|                                                                                          | MicroPoint™                | 100           | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1152        |
| <b>ART® 200</b>                                                                          |                            |               |         |                                   |           |                 |
|        | Low Retention              | 200           | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1165        |
|                                                                                          | Standard                   | 200           | +       | 10 racks, 96 each (Reload insert) | 960       | 732-1153        |
| <b>ART® 1000XL</b>                                                                       |                            |               |         |                                   |           |                 |
|        | Low Retention              | 1000          | +       | 8 racks, 96 each (Reload insert)  | 768       | 732-1172        |
|                                                                                          | Standard                   | 1000          | +       | 8 racks, 96 each (Reload insert)  | 768       | 732-1154        |
| <b>Description</b>                                                                       |                            |               |         |                                   | <b>Pk</b> | <b>Cat. No.</b> |
| <b>Accessories</b>                                                                       |                            |               |         |                                   |           |                 |
| Empty reload rack, sterile, small, for 10 µl and 20 µl tips (insert not included)        |                            |               |         |                                   | 10        | 732-1171        |
| Empty reload rack, sterile, medium, for 20 µl, 200 and 300 µl tips (insert not included) |                            |               |         |                                   | 10        | 732-1170        |
| Empty reload rack, sterile, large, for 1000 µl tips (insert not included)                |                            |               |         |                                   | 8         | 732-1169        |



## VWR COLLECTION CATALOGUE

10 000 equipment, instrument and consumable products

Request your copy from your local VWR sales office or [vwr.com](http://vwr.com)