

Continued from previous page

Capacity (ml)	Dosage	Cap	Packed	Neck Ø ext. (mm)	WxDxH (mm)	Pk	Cat. No.
Rectangular bottles, wide neck							
1000	Sodium thiosulphate (20 mg/l)	Tamper evident screw cap with wad	Carton of 85	55	92x77x187	85	331-0052
1000	Sodium thiosulphate (120 mg/l)	Tamper evident screw cap with wad	Carton of 85	55	92x77x187	85	331-0053
2000	Sodium thiosulphate (20 mg/l)	Tamper evident screw cap with wad	Carton of 50	58	110x93x237	50	331-0054



Bottles, narrow mouth, with screw cap



PET, clear or green

For basic sample collecting. Non sterile, lightweight, round.

- 28 mm neck for all bottles
- Green coloured to protect light sensitive material
- Manufactured from food grade PET

Capacity (ml)	Colour	Pk	Cat. No.
250	Clear	200	331-0094
500	Clear	100	331-0095
1000	Clear	60	331-0096
1000	Green	60	331-0099
2000	Clear	20	331-0098



Bottles, wide mouth, with screw cap



PET, clear

For basic sample collecting. Non sterile, lightweight, round or square.

- 34 mm neck for all bottles
- Tamper evident screw caps
- Manufactured by food grade PET

Capacity (ml)	Type	Cap colour	Form	Pk	Cat. No.
250	Tamper evident screw cap	Blue	Square	200	331-0161
500	Tamper evident screw cap	Blue	Square	100	331-0155
1000	Tamper evident screw cap	Blue	Square	70	331-0167
250	Tamper evident screw cap	Orange	Round	200	331-0165
500	Tamper evident screw cap	Orange	Round	100	331-0166

ALL YOU NEED FOR FOOD ANALYSIS

Specialist publication for all those working in food science

Request your copy from your local VWR sales office or vwr.com