



## Personal eyewash Cederroth



Sterile, buffered, isotonic sodium chloride solution for rinsing eyes following an accident. Available as a single use, pocket version (235 ml) or wall-mounted (500 ml). The bottle opens automatically when twisted out of the wall bracket, and must therefore be used immediately.

- Quick and easy to open
- Eye cup protected from dust
- Generous flow of fluid for effective rinsing - the eye cup on the bottle channels the fluid to the eye, avoiding unnecessary spills
- Neutralises the pH value
- Expiry date: Four years if container not opened

| Description   | Pk | Cat. No. |
|---|----|----------|
| Eyewash solution, sterile, pocket model, 235 ml   | 1  | 111-0164 |
| Eyewash solution, sterile, 500 ml, twin pack  | 2  | 118-0153 |
| Wall bracket for 500 ml bottle*   | 1  | 118-0001 |
| Strap for wall bracket, for extra support for the bottle when located in vibrating environments | 1  | 118-0248 |
| Holster for pocket model, green, with belt loop, 8x4,5x16 mm                                    | 1  | 111-0165 |

\* Bottle holder is supplied without eyewash

## Eyewash station Cederroth



Wall-mountable eyewash station with two 500 ml bottles which open automatically when they are twisted out of their holders. The solution works to neutralise eye pH, i.e. it has a neutralising effect on any alkali or acid splashes, giving a quicker result than a normal sodium chloride solution. Rinsing time: ~1,5 minutes.

- Contains Salvequick plaster dispenser (incl. 40 sterile fabric and 45 sterile plastic plasters), eyewash instructions and a key
- Features extra space for other accessories, e.g. Bloodstopper dressings, protective glasses or gloves
- Transparent door protects the products against dust and dirt and makes it easy to see when a refill is needed

WxDxH: 290x120x560 mm  
Shelf life: 4,5 years

**Ordering information:** The eyewash bottles are included. The Salvequick plaster dispenser is included as well but can also be ordered separately for refilling.

| Description   | Pk | Cat. No. |
|---|----|----------|
| Eyewash station with two 500 ml bottles and plaster dispenser | 1  | 118-0277 |

  

| Description  | Pk | Cat. No. |
|--|----|----------|
| <b>Accessories</b>   |    |          |
| Salvequick plaster dispenser (incl. 40 sterile fabric and 45 sterile plastic plasters) | 1  | 118-0159 |
| Eyewash solution, sterile, 500 ml, twin pack   | 2  | 118-0153 |



## Eye rinse Cederroth



This eyewash solution contains physiological sodium chloride solution (NaCl 0,9%) for rinsing dust and dirt out of the eyes.

- Sterile, for single use
- Relieves irritation in dry eyes
- Vial is easily opened and ready to use by twisting off the top

**Delivery information:** Supplied in a practical cardboard dispenser (incl. 10x 20 ml vials) for wall-mounting, or in a box (incl. 12x 20 ml vials) which fits in the first aid station.

| Description  | Pk | Cat. No. |
|--|----|----------|
| Eyewash solution, 10 vials of 20 ml, packed in dispenser (25x120x202 mm) | 1  | 118-0243 |
| Eyewash solution, 12 vials of 20 ml, packed in box (70x63x118 mm)        | 12 | 118-0244 |