

Beaker and flask holders

Bel-Art Products



Lead ring, with red vinyl coating and four stainless steel brackets

Holds beakers and flasks firmly in a fixed position. Ideal for use in water baths as it allows liquid to circulate beneath the vessel base.

- Good chemical resistance
- Good thermal resistance
- Easy, simple handling

For	I-Ø bottom (mm)	Weight (g)	Pk	Cat. No.
100 ml beaker	60	208	1	214-1935
250 ml beaker	80	340	1	214-1936

Ring weights for vessels

Heathrow Scientific



Lead core, with blue vinyl coating

Open or closed design. Open design, also suitable for vessels with lateral connections. For stabilising glass or plastic vessels on shakers or in water and oil baths, etc. The weights prevent tipping. The coating prevents the vessels from being scratched or causing damage to other surfaces.

- Good chemical resistance
- High temperature resistance
- Available in various sizes for use with a range of vessel sizes

For (ml)	Ø ext. (mm)	Ø int. (mm)	Weight (g)	Pk	Cat. No.
Closed					
Flasks 125 - 500	75	48	270	1	214-1942
Flasks 250 - 1000	90	51	500	1	214-1943
Flasks 500 - 2000	102	57	700	1	214-1944
Flasks 1000 - 4000	115	70	840	1	214-1945
Open					
Flasks 125 - 500	75	42	200	1	214-1891
Flasks 250 - 1000	90	51	480	1	214-1892
Flasks 500 - 2000	102	54	660	1	214-1893
Flasks 1000 - 4000	115	66	930	1	214-1894

Ring weights for vessels

Bel-Art Products



Lead core, with red Vikem® vinyl coating

- Weight prevents tipping and floating in water baths
- Non slip Vikem® coating prevents damage to surfaces
- Good chemical resistance

For (ml)	Ø ext. (mm)	Ø int. (mm)	Weight (g)	Pk	Cat. No.
Cylinders 5 - 10	51	19	165	1	214-1941