

Continued from previous page

Description	Range (°C)	Accuracy (°C)	Division (°C)	Length (mm)	Pk	Cat. No.
For freezers	-25...-5	±0,2	0,1	210	1	620-0837
For refrigerators	-2...+10	±0,5	0,1	180	1	620-0838
For incubators	+15...+30	±0,2	0,1	180	1	620-0839
For incubators	+25...+45	±0,2	0,1	190	1	620-0840
For ovens	+95...+115	±0,2	0,1	200	1	620-0841



## FRIO-Temp® precision thermometers



### Enviro-Safe® liquid or red spirit

These glass thermometers can be used to take readings in freezers, refrigerators, incubators, ovens and blood banks.

- Enviro-Safe® non toxic and biodegradable liquid
- Individual serial number
- Wide choice

Traceable to NIST/DKD

**Delivery information:** Supplied with a statement of accuracy plus Velcro® and magnetic attachments in recyclable packaging.

Description	Range (°C)	Accuracy (°C)	Division (°C)	Length (mm)	Pk	Cat. No.
<b>Enviro-Safe®</b>						
For freezers	-30...0	±0,5	0,5	135	1	620-0854
For refrigerators	-5...+15	±0,5	0,5	125	1	620-0902
For blood banks	-5...+20	±0,5	0,5	135	1	620-0848
For incubators	+15...+50	±0,5	0,5	135	1	620-0877
For incubators	+15...+50	±0,5	0,5	135	1	620-0878
For incubators	+15...+50	±0,5	0,5	135	1	620-0800
For ovens	+20...+130	±1,0	1,0	145	1	620-0850
For ovens	+50...+200	±2,0	2,0*	145	1	620-0852
<b>Red spirit</b>						
For ultra-low freezers	-90...+20	±1,0	1,0**	145	1	620-0855
For freezers	-30...+1	±0,5	0,5	135	1	620-0843
For refrigerators	-5...+15	±0,5	0,5	125	1	620-0845
For blood banks	-5...+20	±0,5	0,5	135	1	620-0842
For incubators	+15...+50	±0,5	0,5	135	1	620-0856
For ovens	+20...+130	±1,0	1,0	145	1	620-0844
For ovens	+50...+220	±2,0	2,0*	145	1	620-0886

\* 3,0 °C (>130 °C) \*\* 2,0 °C (<-50 °C)



## FRIO-Temp® precision thermometers, VWR



### Black spirit

These thermometers are designed for applications in freezers, refrigerators, incubators, ovens and blood banks. To use, simply screw thermometer with cap into bottle and attach to test sample.

- Very accurate
- No hazard of mercury
- Individual serial number

Traceable to standard NIST/DKD.

**Delivery information:** Supplied with a statement of accuracy plus Velcro® and magnetic attachments in recyclable packaging.