

GC injection port septa



Long-life GC-septa for inlet temperature of up to 340 °C.

- Low bleeding
- Easy exchange, no sticking to hot surfaces
- Good penetration and sealing properties

Description	Ø (mm)	Pk	Cat. No.
Shimadzu plug, transparent	-	25	548-0181
Universal high temperature injection port septa, transparent	9,5	25	548-0182
Universal high temperature injection port septa, transparent	10	25	548-0183
Universal high temperature injection port septa, transparent	11	25	548-0184
Universal high temperature injection port septa, transparent	12,5	25	548-0185
Universal high temperature injection port septa, transparent	17	25	548-0186

Gas-clean filters Agilent Technologies



Clean gases reduce the risk of column damage and improve accuracy and sensitivity. These unbreakable and self-indicating filters give high sensitivity to oxygen and moisture and can be installed in seconds without tools.

Description	Pk	Cat. No.
Oxygen filter	1	558-0235
Moisture filter	1	569-7971
Charcoal filter	1	558-0236
GC/MS filter	1	558-0237
Installation kit GC/MS	1	558-0228
Connecting unit for 4 filters (1/8" tube)	1	558-0229
Connecting unit for 2 filters (1/4" tube)	1	558-0049
Connecting unit for 2 filters (1/8" tube)	1	558-0231
Connecting unit for 1 filter (1/4" tube)	1	558-0233
Connecting unit for 1 filter (1/8" tube)	1	558-0048
Connecting unit for 4 filters (1/4" tube)	1	558-0051
FID-kit 1/8" includes 2 charcoal, 1 oxygen, 1 moisture filter, connecting unit for 4 filters	1	558-0230
TCD-kit 1/8" includes 1 oxygen, 1 moisture filter, connecting unit for 2 filters	1	558-0232
Wall mounting bracket/single connecting unit	1	558-0053
FID-kit 1/4" includes 2 charcoal, 1 oxygen, 1 moisture filter, connecting unit for 4 filters	1	558-0234

Capillary columns, BPX70 SGE (Trajan)



70% cyanopropyl polysilphenylene-siloxane

Polarity has been designed for fatty acid methyl ester (FAME) analysis.

Max. temperature: 260 °C

Suitable replacement for: CP-Sil88, DB-23, HP-23, Rtx-2330, SP-2330, SP-2380.

Length (m)	Ø int. (mm)	Film (µm)	Pk	Cat. No.
10	0,1	0,2	1	568-0170
15	0,25	0,25	1	568-0180
25	0,22	0,25	1	558-0332
25	0,32	0,25	1	558-0334
30	0,22	0,25	1	568-0174
30	0,25	0,25	1	558-0252
30	0,32	0,25	1	568-0176
30	0,53	0,5	1	568-0179
50	0,22	0,25	1	558-0333
50	0,32	0,25	1	568-0172
60	0,22	0,25	1	568-0175
60	0,25	0,25	1	558-0335
60	0,32	0,25	1	568-0177
120	0,25	0,25	1	568-0181

Further chromatography consumables, instruments and accessories are available in our AYN for Chromatography catalogue.