

Dispenser tips, Combitips® Advanced
Eppendorf



Combitips® Advanced



Rack for Combitips®

For use with Multipette®.

- Ergonomic design with unique funnel shape prevents damage to gloves and ensures comfortable handling
- Individually colour coded for quick identification
- Elongated tips (for 2,5, 5 and 10 ml) allow complete emptying of all common tubes preventing sample loss

Combitips® advanced and Combitips® advanced PCR Clean are supplied in a dispenser box (100 Combitips® in 4 bags, PCR Clean with zip-lock bag); 25/50 ml Combitips® come in 25 pcs. per box (each containing an adapter).

Combitips® Biopur® are individually packed in single blister and additional plastic bag for cleanroom use.

Capacity (ml)	Sterile	Packed	Pk	Cat. No.
Combitips® Advanced				
0,1	-	4 bags, 25 each	100	613-2058
0,2	-	4 bags, 25 each	100	613-2059
0,5	-	4 bags, 25 each	100	613-2060
1,0	-	4 bags, 25 each	100	613-2061
2,5	-	4 bags, 25 each	100	613-2062
5,0	-	4 bags, 25 each	100	613-2063
10	-	4 bags, 25 each	100	613-2064
25	-	4 bags, 25 each	100	613-2065
50	-	4 bags, 25 each	100	613-2066
Combitips® Advanced PCR Clean				
0,1	-	4 bags, 25 each	100	613-2076
0,2	-	4 bags, 25 each	100	613-2077
0,5	-	4 bags, 25 each	100	613-2078
1,0	-	4 bags, 25 each	100	613-2079
2,5	-	4 bags, 25 each	100	613-2080
5,0	-	4 bags, 25 each	100	613-2081
10	-	4 bags, 25 each	100	613-2082
25	-	4 bags, 25 each	100	613-2083
50	-	4 bags, 25 each	100	613-2084
Combitips® Advanced Biopur®				
0,1	+	Individual blister	100	613-2067
0,2	+	Individual blister	100	613-2068
0,5	+	Individual blister	100	613-2069
1,0	+	Individual blister	100	613-2070
2,5	+	Individual blister	100	613-2071
5,0	+	Individual blister	100	613-2072
10	+	Individual blister	100	613-2073
25	+	Individual blister	100	613-2074
50	+	Individual blister	100	613-2075
Description			Pk	Cat. No.
Accessories				
Adapter 25 ml for Combitips® Advanced			1	613-2085
Adapter 50 ml for Combitips® Advanced			1	613-2086
Adapter 25 ml for Combitips® Advanced Biopur®			7	613-2087
Adapter 50 ml for Combitips® Advanced Biopur®			7	613-2088
Rack for 8 Combitips® Advanced			1	613-2089



Pipette tips, Finntip® PDP
Thermo Scientific



Autoclavable

- Capillary unit and positive displacement piston
- Recommended for PCR, DNA sequencing and all delicate pipetting tasks (viscous or foaming liquids, etc.)
- For use with Finnpiettes® PDP

Capacity (µl)	Sterile	Packed	Pk	Cat. No.
0,5 - 25	-	10 racks, 96 each	960	613-2630
20 - 200	-	10 racks, 96 each	960	613-2631
0,5 - 25	-	Bulk	960	613-2632
20 - 200	-	Bulk	960	613-2633