► Thermo Scientific Sterilin Swabs, Plain





Sterilin Plain Swabs are supplied in either a peel pouch or in a tube for ease of handling and transportation.

details

- Primary sample collection device for bacteria
- · Color-coded caps for ease of identification

Ordering Information: Available in UK only

Certifications: CE marked as Class IIa in accordance with the Medical Device Directive 93/42/EEC (for transient invasive use)

Notes: Must be processed as soon as possible after sampling (ideally within two-hour period) to avoid dehydration of the sample and possible non-recovery of pathogens

Sterilin Swabs, Plain

Applicator/Tip	Cap Color	Use	Sterility	Inner Pack Qty.	Case Qty.
ibe					
Wood/Cotton	Red	General purpose	IRR	1	100
Plastic/Rayon	White	General purpose	IRR	1	100
Aluminum/Rayon	Orange	Ear, nose, vaginal or urethral	IRR	1	100
Twisted wire/Rayon	Blue	Pernasal and nasopharyngeal	IRR	1	100
Plastic/Cotton, Vegetal protein coated	Yellow	Vegetal protein coating reduces overgrowth	IRR	1	100
Wood/Cotton, Charcoal coated	Black	Charcoal neutralizes inhibitory effects of cotton	IRR	1	100
eel pouch					
Wood/Cotton	N/A	General purpose	E0	1	1,000
Wood/Cotton	N/A	General purpose	IRR	100	2,500
	Wood/Cotton Plastic/Rayon Aluminum/Rayon Twisted wire/Rayon Plastic/Cotton, Vegetal protein coated Wood/Cotton, Charcoal coated eel pouch Wood/Cotton	Wood/Cotton Red Plastic/Rayon White Aluminum/Rayon Orange Twisted wire/Rayon Blue Plastic/Cotton, Vegetal protein coated Yellow Wood/Cotton, Charcoal coated Black eel pouch Wood/Cotton N/A	Wood/Cotton Red General purpose Plastic/Rayon White General purpose Aluminum/Rayon Orange Ear, nose, vaginal or urethral Twisted wire/Rayon Blue Pernasal and nasopharyngeal Plastic/Cotton, Vegetal protein coated Yellow Vegetal protein coating reduces overgrowth Wood/Cotton, Charcoal coated Black Charcoal neutralizes inhibitory effects of cotton eel pouch Wood/Cotton N/A General purpose	Wood/Cotton Red General purpose IRR Plastic/Rayon White General purpose IRR Aluminum/Rayon Orange Ear, nose, vaginal or urethral IRR Twisted wire/Rayon Blue Pernasal and nasopharyngeal IRR Plastic/Cotton, Vegetal protein coated Yellow Vegetal protein coating reduces overgrowth IRR Wood/Cotton, Charcoal coated Black Charcoal neutralizes inhibitory effects of cotton IRR wood/Cotton N/A General purpose EO	Wood/Cotton Red General purpose IRR 1 Plastic/Rayon White General purpose IRR 1 Aluminum/Rayon Orange Ear, nose, vaginal or urethral IRR 1 Twisted wire/Rayon Blue Pernasal and nasopharyngeal IRR 1 Plastic/Cotton, Vegetal protein coated Yellow Vegetal protein coating reduces overgrowth IRR 1 Wood/Cotton, Charcoal coated Black Charcoal neutralizes inhibitory effects of cotton IRR 1 wood/Cotton N/A General purpose EO 1