PCR plates, 96-well, ABgene® SuperPlate™ Thermo Scientific




PF

The ABgene® SuperPlate™ provides four times more rigidity for superior robotic handling. Three versions are available, all with ultra-rigid SBS footprint for compatibility with automated systems, and thin walled round bottom wells for optimal heat transfer and maximum sample recovery.

- Enhanced thermal stability to minimise warping during thermal cycling and heat sealing
- Black lettering for quick referencing of wells
- Maximum well volume 0,2 ml when used with adhesive heat seals
- Skirted version is stackable and compatible with automated systems
- Available in white for qPCR applications

| Description | Colour | Pk | Cat. No. |
|-------------------|---------|----|----------|
| Raised semi-skirt | Blue | 25 | 732-1057 |
| Raised semi-skirt | Purple | 25 | 732-1058 |
| Raised semi-skirt | White | 25 | 731-0305 |
| Raised semi-skirt | Yellow | 25 | 732-1059 |
| Raised semi-skirt | Natural | 25 | 732-1056 |
| Semi-skirted | Blue | 25 | 732-1061 |
| Semi-skirted | Purple | 25 | 732-1062 |
| Semi-skirted | White | 25 | 732-1063 |
| Semi-skirted | Yellow | 25 | 732-1064 |
| Semi-skirted | Natural | 25 | 732-1060 |
| Skirted | Blue | 25 | 732-1050 |
| Skirted | Purple | 25 | 732-1051 |
| Skirted | White | 25 | 732-1052 |
| Skirted | Yellow | 25 | 732-1053 |
| Skirted | Natural | 25 | 732-1049 |

PCR plates, 96-well, fully skirted, low profile, ABgene® Thermo-Fast® Thermo Scientific




PP

The low profile of the plate decreases 'dead space' between the heated lid of the thermal cycler and the sample. This eliminates condensation on the side-wall of the tube, preventing reduction in PCR volume and increasing the efficiency of the reaction.

- Suitable for use with automated systems
- Maximum well volume 0,2 ml when used with adhesive and heat seals
- Eight holes in skirt aid plate positioning and removal from thermal cycler block
- Clear well bottom for sample visibility
- Certified free from DNase, RNase and human genomic DNA

| Description | Colour | Pk | Cat. No. |
|--|---------|----|----------|
| PCR plate, 96-well, fully skirted, low profile | White | 25 | 732-0391 |
| PCR plate, 96-well, fully skirted, low profile | Natural | 25 | 732-4888 |

PCR plates, 96-well, rigid, semi-skirted, ABgene® Thermo-Fast® Thermo Scientific



PΡ

A rigid version of the standard semi-skirted plate for applications where extra plate rigidity is required. Colour: Clear

| Description | Pk | Cat. No. |
|---|----|----------|
| PCR plate, 96-well, rigid, semi-skirted | 25 | 732-4924 |